

Paths of analysis*

L11_DIA

Synthia

October 11, 2022

1 Analysis parameters

Analysis type: Automatic Retrosynthesis

Rules: none selected

Filters: Tunnels, FGI, FGI with protections

Max. paths returned: 50

Max. iterations: 2000

Commercial:

1. Max. molecular weight - 1000 g/mol
2. Max. price - 1500 \$/g

Published:

1. Max. molecular weight - 1000 g/mol
2. Popularity - 5

My Stockroom:

1. Max. molecular weight - 1000 g/mol

Reaction scoring formula: $\text{TUNNEL_COEF} * \text{FGI_COEF} * \text{STEP} * 20 + 1000000 * (\text{CONFLICT} + \text{NON_SELECTIVITY} + \text{FILTERS} + \text{PROTECT})$

Chemical scoring formula: $\text{SMALLER}^3, \text{SMALLER}^{1.5}$

Min. search width: 400

Max. reactions per product: 60

Strategies: none selected

*The results stated herein were generated using the proprietary platform owned and maintained by Grzybowski Scientific Inventions, Inc., a subsidiary of Merck KGaA, Darmstadt Germany. The results are provided on an as is basis, and shall be used solely in connection with the rights afforded in the license agreement and for no other purpose.

JSON Parameters: $\{\}$

2 paths found. *Paths are sorted by score. Reactions are sorted in appearance order for each path.*

Score: 360.85




Figure 1: Outline of path 1

1. Methyl chloroformate - *available at Sigma-Aldrich*
2. 3-bromo-2-methoxy-propene


1. methyl 3-methoxy-3-butenolate

Protections: none

Reference: [10.1016/0040-4039\(94\)85361-4](#) and [10.1016/0040-4039\(96\)00258-4](#)
and [10.1021/jo3005556](#) AND [10.1016/0040-4039\(96\)00689-2](#)

Retrosynthesis ID: 2530

2.1.2 Oxidation of primary alcohols with DMP



Substrates:

1. 2,3-Dibromoallyl alcohol - *available at Sigma-Aldrich*

Products:

1. 2,3-dibromo-propenal


Typical conditions: DMP.DCM.0-25 C

Protections: none

Reference: [10.1016/j.bmc.2020.115469](#) p. 3, 9 and [10.1021/acs.jmedchem.8b01878](#) SI p. S43

Retrosynthesis ID: 50426

2.1.3 Condensation of esters with aldehydes/ketones



Substrates:

1. methyl 3-methoxy-3-butenate
2. 2,3-dibromo-propenal

Products:

1. C=C(OC)C(=CC(Br)=CBr)C(=O)OC


Typical conditions: LDA.THF

Protections: none

Reference: [10.1021/op040006z](#) AND [10.1016/j.bmcl.2005.10.104](#) AND

Retrosynthesis ID: 14983

2.1.4 Synthesis of ketones and aldehydes from enol ethers



Substrates:

1. C=C(OC)C(=CC(Br)=CBr)C(=O)OC

Products:

1. COC(=O)C(=CC(Br)=CBr)C(C)=O


Typical conditions: [H+].THF

Protections: none

Reference: [10.1081/SCC-120023437](#) AND [10.1016/j.bmcl.2007.11.020](#) AND [10.1016/j.tet.2011.03.084](#) AND [10.1021/ja00270a023](#) AND [10.1055/s-1994-25424](#) AND

Retrosynthesis ID: 14844

2.1.5 Shapiro reaction followed by DMF addition



Substrates:

1. 3-(p-methoxybenzyloxy)propanal
2. Dimethylformamide - *available at Sigma-Aldrich*

Products:

1. COc1ccc(COCC=CC=O)cc1


Typical conditions: 1. TsNH₂NH₂. 2. Mes₂Mg.LiCl. THF. heating then DMF. 3. NaBH₄

Protections: none

Reference: [10.1021/jo901926z](#) and [10.1002/chem.201303586](#) (SI p.7) and [10.1021/jo015699l](#)

Retrosynthesis ID: 9990428

2.1.6 Deprotection of PMB ethers



Substrates:

1. COc1ccc(COCC=CC=O)cc1

Products:

1. 4-hydroxy-but-2-enal


Typical conditions: DDQ.H₂O.DCM or CAN.H₂O.ACN

Protections: none

Reference: [10.1080/00397911.2019.1639757](#) and [10.1039/C9OB01504C](#) and [10.1080/00397911.2019.1660900](#) and [10.1039/C9OB00556K](#)

Retrosynthesis ID: 31010150

2.1.7 Acid catalyzed transesterification



Substrates:

1. 4-hydroxy-but-2-enal
2. COC(=O)C(=CC(Br)=CBr)C(C)=O

Products:

1. CC(=O)C(=CC(Br)=CBr)C(=O)OCC=CC=O


Typical conditions: H⁺

Protections: none

Reference: [10.1021/cr00020a004](#)

Retrosynthesis ID: 50438

2.1.8 Diels-Alder



Substrates:

1. CC(=O)C(=CC(Br)=CBr)C(=O)OCC=CC=O

Products:

1. CC(=O)C12C=C(Br)C(Br)C(C=O)C1COC2=O


Typical conditions: Lewis acid or chiral Lewis acid. Solvent.

Protections: none

Reference: DOI: [10.1002/1521-3773\(20020517\)41:10<1668::AID-ANIE1668>3.0.CO;2-Z](#) AND [10.1021/ja062508t](#)

Retrosynthesis ID: 18116

2.1.9 Takai olefination



Substrates:

- 1,1-Diiodoethane - *available at Sigma-Aldrich*
- CC(=O)C12C=C(Br)C(Br)C(C=O)C1COC2=O

Products:

- C/C=C/C1C(Br)C(Br)=CC2(C(C)=O)C(=O)OCC12


Typical conditions: CrCl2.THF.DMF

Protections: none

Reference: [10.1021/ja00283a046](#) and [10.1021/ja00237a081](#)

Retrosynthesis ID: 10942

2.1.10 Alkylation of amines with alkyl bromides



Substrates:

- 2-(p-Toluenesulfonylamino)phenylboronic acid pinacol ester - *available at Sigma-Aldrich*
- C/C=C/C1C(Br)C(Br)=CC2(C(C)=O)C(=O)OCC12

Products:

1. C/C=C/C1C(N(c2ccccc2B2OC(C)(C)C(C)(C)O2)S(=O)(=O)c2ccc(C)cc2)C(Br)=CC2(C(C)=O)C(=O)C2


Typical conditions: K₂CO₃ or other base

Protections: none

Reference: [10.1016/j.tetlet.2007.09.110](https://doi.org/10.1016/j.tetlet.2007.09.110)

Retrosynthesis ID: 7668

2.1.11 Suzuki coupling of arylboronic pinacol esters with vinyl Bromides



Substrates:

1. C/C=C/C1C(N(c2ccccc2B2OC(C)(C)C(C)(C)O2)S(=O)(=O)c2ccc(C)cc2)C(Br)=CC2(C(C)=O)C(=O)C2

Products:

1. C/C=C/C1C2C(=CC3(C(C)=O)C(=O)OCC13)c1ccccc1N2S(=O)(=O)c1ccc(C)cc1

Typical conditions: Pd catalyst.base.solvent

Protections: none

Reference: [10.1021/cr00039a007](https://doi.org/10.1021/cr00039a007) and [10.1007/3418_2012_32](https://doi.org/10.1007/3418_2012_32) and [10.1021/cr0505268](https://doi.org/10.1021/cr0505268) and [10.1016/j.jfluchem.2016.01.018](https://doi.org/10.1016/j.jfluchem.2016.01.018) and [10.1039/C3CS60197H](https://doi.org/10.1039/C3CS60197H)

Retrosynthesis ID: 10695

2.2 Path 2

Score: 369.43





Figure 2: Outline of path 2

2.2.1 Bromination of Alkynes



Substrates:

1. 3-Bromo-1-propyne - *available at Sigma-Aldrich*

Products:

1. 1,2,3-tribrom-propen


Typical conditions: NBS. THF. 80C

Protections: none

Reference: DOI: [10.1016/j.tetlet.2011.06.047](https://doi.org/10.1016/j.tetlet.2011.06.047) or DOI: [10.1055/s-2006-941558](https://doi.org/10.1055/s-2006-941558) or DOI: [10.1021/jo011016q](https://doi.org/10.1021/jo011016q)

Retrosynthesis ID: 8354

2.2.2 Aldehyde Formation



Substrates:

1. methyl 4-hydroxybut-2-ynoate - *available at Sigma-Aldrich*

Products:

1. 4-hydroxy-but-2-ynal


Typical conditions: DIBAL.solvent e.g. DCM

Protections: none

Reference: [10.1039/C39940000483](#) and [10.1039/C3CC47867J](#) and [10.1021/jo00222a054](#) and [10.1021/ja9934908](#) and [10.1021/jo902426z](#)

Retrosynthesis ID: 28551

2.2.3 HWE/Wittig Olefination



Substrates:

1. 1,2,3-tribrom-propen
2. Methyl 2-ketobutyrate - *available at Sigma-Aldrich*

Products:

1. CCC(=CC(Br)=CBr)C(=O)OC


Typical conditions: 1.PPh₃ or trialkylphosphite.2.base.aldehyde


Protections: none

Reference: [10.1002/anie.200705005](#) and [10.1021/ol052106a](#) and [10.1021/jo00075a064](#) and [10.1021/ol3027297](#)

Retrosynthesis ID: 24425

2.2.4 Allylic Oxidation of Alkenes





Substrates:

1. CCC(=CC(Br)=CBr)C(=O)OC

Products:

1. COC(=O)C(=CC(Br)=CBr)C(C)=O


Typical conditions: tBuOOH.Pd(OH)₂/C or PhI(OAc)₂ or SeO₂

Protections: none

Reference: [10.1021/ja0340735](#) and [10.1021/ol100603q](#) and [10.1016/j.tetlet.2016.05.063](#) (Scheme 2)

Retrosynthesis ID: 2583

2.2.5 Acid catalyzed transesterification



Substrates:

1. 4-hydroxy-but-2-ynal
2. COC(=O)C(=CC(Br)=CBr)C(C)=O

Products:

1. CC(=O)C(=CC(Br)=CBr)C(=O)OCC#CC=O


Typical conditions: H⁺

Protections: none

Reference: [10.1021/cr00020a004](#)

Retrosynthesis ID: 50438

2.2.6 Reduction of alkynes to alkenes



Substrates:

1. CC(=O)C(=CC(Br)=CBr)C(=O)OCC#CC=O

Products:

1. CC(=O)C(=CC(Br)=CBr)C(=O)OCC=CC=O


Typical conditions: H₂.Lindlar's catalyst

Protections: none

Reference: [10.1021/ja054487t](#) (SI) AND [10.1021/jm9810912](#)

Retrosynthesis ID: 14627

2.2.7 Diels-Alder



Substrates:

1. CC(=O)C(=CC(Br)=CBr)C(=O)OCC=CC=O

Products:

1. CC(=O)C12C=C(Br)C(Br)C(C=O)C1COC2=O


Typical conditions: Lewis acid or chiral Lewis acid. Solvent.

Protections: none

Reference: DOI: [10.1002/1521-3773\(20020517\)41:10<1668::AID-ANIE1668>3.0.CO;2-Z](https://doi.org/10.1002/1521-3773(20020517)41:10<1668::AID-ANIE1668>3.0.CO;2-Z) AND [10.1021/ja062508t](https://doi.org/10.1021/ja062508t)

Retrosynthesis ID: 18116

2.2.8 Takai olefination



Substrates:

- 1,1-Diiodoethane - *available at Sigma-Aldrich*
- CC(=O)C12C=C(Br)C(Br)C(C=O)C1COC2=O

Products:

- C/C=C/C1C(Br)C(Br)=CC2(C(C)=O)C(=O)OCC12


Typical conditions: CrCl₂.THF.DMF

Protections: none

Reference: [10.1021/ja00283a046](https://doi.org/10.1021/ja00283a046) and [10.1021/ja00237a081](https://doi.org/10.1021/ja00237a081)

Retrosynthesis ID: 10942

2.2.9 Alkylation of amines with alkyl bromides




Substrates:

1. 2-(p-Toluenesulfonylamino)phenylboronic acid pinacol ester - *available at Sigma-Aldrich*
2. C/C=C/C1C(Br)C(Br)=CC2(C(C)=O)C(=O)OCC12

Products:

1. C/C=C/C1C(N(c2ccccc2B2OC(C)(C)C(C)(C)O2)S(=O)(=O)c2ccc(C)cc2)C(Br)=CC2(C(C)=O)C(=O)OCC12

Typical conditions: K₂CO₃ or other base**Protections:** none**Reference:** [10.1016/j.tetlet.2007.09.110](https://doi.org/10.1016/j.tetlet.2007.09.110)**Retrosynthesis ID:** 7668**2.2.10 Suzuki coupling of arylboronic pinacol esters with vinyl Bromides****Substrates:**

1. C/C=C/C1C(N(c2ccccc2B2OC(C)(C)C(C)(C)O2)S(=O)(=O)c2ccc(C)cc2)C(Br)=CC2(C(C)=O)C(=O)OCC12

Products:

1. C/C=C/C1C2C(=CC3(C(C)=O)C(=O)OCC13)c1ccccc1N2S(=O)(=O)c1ccc(C)cc1

Typical conditions: Pd catalyst.base.solvent**Protections:** none**Reference:** [10.1021/cr00039a007](https://doi.org/10.1021/cr00039a007) and [10.1007/3418_2012_32](https://doi.org/10.1007/3418_2012_32) and [10.1021/cr0505268](https://doi.org/10.1021/cr0505268) and [10.1016/j.jfluchem.2016.01.018](https://doi.org/10.1016/j.jfluchem.2016.01.018) and [10.1039/C3CS60197H](https://doi.org/10.1039/C3CS60197H)

Retrosynthesis ID: 10695