Paths of analysis*

Synthia

October 10, 2022

1 Analysis parameters

Analysis type: Automatic Retrosynthesis

Rules: none selected

Filters: Exclude Diastereoselecitve reactions, Tunnels, FGI, FGI with protec-

tions

Max. paths returned: 50

Max. iterations: 2000

Commercial:

1. Max. molecular weight - 1000 g/mol

2. Max. price - 1500 \$/g

Published:

- 1. Max. molecular weight 1000 g/mol
- 2. Popularity 5

My Stockroom:

1. Max. molecular weight - 1000 g/mol

Reaction scoring formula: TUNNEL_COEF*FGI_COEF*STEP*20+1000 000*(CONFLICT+NON SELECTIVITY+FILTERS+PROTECT)

Chemical scoring formula: SMALLER^ 3,SMALLER^ 1.5

Min. search width: 400

Max. reactions per product: 60

^{*}The results stated herein were generated using the proprietary platform owned and maintained by Grzybowski Scientific Inventions, Inc., a subsidiary of Merck KGaA, Darmstadt Germany. The results are provided on an as is basis, and shall be used solely in connection with the rights afforded in the license agreement and for no other purpose.

 ${f Strategies:}\ {f none}\ {f selected}$

FGI Coeff: 0

Tunnels Coeff: 0

JSON Parameters: {}

2 Paths

5 paths found. Paths are sorted by score. Reactions are sorted in appearance order for each path.

2.1 Path 1

Score: 20.00





Figure 1: Outline of path 1

2.1.1 Takai olefination



Substrates:

- 1. 3-Formylphenyl acetate available at Sigma-Aldrich
- 2. Iodoform available at Sigma-Aldrich

Products:

 $1. \ \mathrm{CC(=O)Oc1cccc}(/\mathrm{C=C/I})c1$


 $\textbf{Typical conditions:} \ \mathrm{CrCl2.THF}$

Protections: none

Reference: 10.1021/ja00283a046 and 10.1021/ja00237a081

Retrosynthesis ID: 10497

2.1.2 Hydrolysis of acetates



Substrates:

 $1. \ \mathrm{CC(=O)Oc1cccc}(/\mathrm{C=C/I})c1$

Products:

 $1. \ Oc1cccc(/C{=}C/I)c1$

 $\textbf{Typical conditions:} \ \mathrm{KOH.MeOH}$


Protections: none

Reference: 10.3762/bjoc.10.40 and 10.1016/j.bmc.2009.11.035 and

10.1016/S0040-4020(02)01584-3

Retrosynthesis ID: 32805

2.1.3 Suzuki coupling of alkyl-9-BBNs with vinyl iodides



 $1. \ Oc1cccc(/C{=}C/I)c1$

2. 9-(3,3-pentamethyleneallyl)-9-borabicyclo3.3.1nonane

Products:

 $1. \ \, Oc1cccc(/C=C/CC=C2CCCC2)c1$

 ${\bf Typical\ conditions:}\ {\bf Pd\ catalyst.base.solvent}$

Protections: none

Reference: 10.1021/jo015995y and 10.1016/j.tetlet.2010.11.139 And 10.1021/ol0600741 and 10.1055/s-2002-32602 and 10.1002/anie.200501760

Retrosynthesis ID: 25168

2.2 Path 2

Score: 25.00





Figure 2: Outline of path 2

2.2.1 Synthesis Of Alkyl Iodides Via Appel Reaction



Substrates:

1. (e)-3-(3-methoxyphenyl)-2-propen-1-ol

Products:

1. COc1cccc(/C=C/CI)c1

Typical conditions: Imidazole.PPh3.I2

Protections: none

Reference: 10.1002/1099-0690(200102)2001:3<493::AID-EJOC493>3.0.CO2-B

(compound 20) and 10.1016/j.tet.2014.09.030

Retrosynthesis ID: 9990040

2.2.2 Palladium catalysed alkylation of vinyl iodides

Substrates:

1. COc1cccc(/C=C/CI)c1

2. iodomethylene cyclohexane

Products:


 $1. \ COc1cccc(/C=C/CC=C2CCCC2)c1$

Typical conditions: [Pd].catalyst

Protections: none

Reference: 10.1016/j.bmcl.2005.12.066 and 10.1021/ol052070m and 10.1021/ol5023195 and 10.1002/anie.200703134 and 10.1016/j.bmcl.2005.09.084 and 10.1021/ol0344873

2.2.3 Demethylation of Phenols



Substrates:

1. COc1cccc(/C=C/CC=C2CCCC2)c1

Products:

1. Oc1cccc(/C=C/CC=C2CCCC2)c1

Typical conditions: BBr3.CH2Cl2

Protections: none

Reference: DOI: 10.1021/ja00105a021 and 10.1021/jm00176a011 and 10.1021/jm970277i and 10.1021/ja0106164 and Patent: US2010/16298, 2010,

 $A1, \; page \; 185$

Retrosynthesis ID: 10011837

2.3 Path 3

Score: 31.25





Figure 3: Outline of path 3

2.3.1 Appel Reaction



Substrates:

1. (e)-3-bromocinnamyl alcohol

Products:

 $1. \ BrC/C{=}C/c1cccc(Br)c1$

Typical conditions: PPh3.CBr4


Protections: none

Reference: 10.1021/ja800574m and 10.1016/j.tet.2012.05.010 and

10.1016/j.tet.2004.09.021 (experimental)

Retrosynthesis ID: 9990037

2.3.2 Shapiro reaction followed by alkyl bromide addition



1. Hexahydrobenzaldehyde - available at Sigma-Aldrich

2. BrC/C=C/c1ccc(Br)c1

Products:

1. Brc1cccc(/C=C/CC=C2CCCC2)c1


Typical conditions: 1.TsNH2NH2.2.Mes2Mg.LiCl.THF.heating then alkyl bromide.cooling

Protections: none

Reference: 10.1016/S0040-4039(00)75263-4 and 10.1021/ol300652k and 10.1021/jo015699l

Retrosynthesis ID: 9990458

2.3.3 Br/Li exchange



Substrates:

1. n-BuLi - available at Sigma-Aldrich

2. Brc1cccc(/C=C/CC=C2CCCC2)c1

Products:


1. [Li]c1cccc(/C=C/CC=C2CCCC2)c1

Typical conditions: nBuLi.or.tBuLi.THF.-78C

Protections: none

Reference: 10.1002/ejoc.201101490 and 10.1016/j.tet.2012.03.058 and 10.1016/j.tetlet.2015.01.032 and 10.1021/ja0541175 and 10.1016/j.tetlet.2016.06.123

${\bf 2.3.4}\quad {\bf Addition\ of\ electrophiles\ to\ lithiated\ arenes/heteroarenes}$



Substrates:

1. [Li]c1cccc(/C=C/CC=C2CCCC2)c1

Products:

 $1. \ \, \text{Oc1cccc}(/\text{C=C/CC=C2CCCC2})\text{c1}$

Typical conditions: B(OMe)3 then H2O2.THF

Protections: none

Reference: 10.1039/C7CC09187G (SI) and 10.1002/ejoc.201701142 and

10.1021/acscatal.6b03380 (SI,p.10) and 10.1002/chem.201702143 (SI)

Retrosynthesis ID: 10019525

2.4 Path 4

Score: 31.25





Figure 4: Outline of path 4

2.4.1 Ozonolysis



Substrates:

1. 3-ethenylcyclohexan-1-one - available at Sigma-Aldrich

Products:

1. 3-oxo-cyclohexancarbaldehyd


 $\textbf{Typical conditions:} \ \ O3. MeOH. CH2Cl2. PPh3 \ or \ Me2S. low \ temperature$

Protections: none

Reference: 10.1016/j.tet.2017.03.039

Retrosynthesis ID: 5074

2.4.2 Synthesis Of Alkyl Iodides Via Appel Reaction



1. 2-cyclohexylideneethan-1-ol - available at Sigma-Aldrich

Products:

1. ICC=C1CCCCC1

Typical conditions: Imidazole.PPh3.I2


Protections: none

Reference: 10.1002/1099-0690(200102)2001:3<493::AID-EJOC493>3.0.CO2-B

(compound 20) and 10.1016/j.tet.2014.09.030

Retrosynthesis ID: 9990040

2.4.3 Takai olefination



Substrates:

1. Iodoform - available at Sigma-Aldrich

2. 3-oxo-cyclohexancarbaldehyd

Products:

1. O=C1CCCC(/C=C/I)C1

Typical conditions: CrCl2.THF

Protections: none

Reference: 10.1021/ja00283a046 and 10.1021/ja00237a081

2.4.4 Palladium catalysed alkylation of vinyl iodides

Substrates:

- $1. \ O{=}C1CCCC(/C{=}C/I)C1$
- 2. ICC=C1CCCCC1

Products:

1. O=C1CCCC(/C=C/CC=C2CCCCC2)C1


Typical conditions: [Pd].catalyst

Protections: none

Reference: 10.1016/j.bmcl.2005.12.066 and 10.1021/ol052070m and 10.1021/ol5023195 and 10.1002/anie.200703134 and 10.1016/j.bmcl.2005.09.084 and 10.1021/ol0344873

Retrosynthesis ID: 25162

2.4.5 Synthesis of enones from ketones



Substrates:

 $1. \ O{=}C1CCCC(/C{=}C/CC{=}C2CCCC2)C1$

Products:

1. O=C1C=CCC(/C=C/CC=C2CCCC2)C1

Typical conditions: TMSOTf then oxidant eg PdCl2.O2


Protections: none

Reference: 10.1002/anie.201512005 and 10.3762/bjoc.10.15 and

10.1021/acs.joc.5b02550 and 10.1055/s-2002-32603

Retrosynthesis ID: 10004759

2.4.6 DDQ mediated aromatization



Substrates:

 $1. \ O{=}C1C{=}CCC(/C{=}C/CC{=}C2CCCC2)C1$

Products:

 $1. \ \, Oc1cccc(/C{=}C/CC{=}C2CCCC2)c1$

Typical conditions: DDQ

Protections: none

Reference: 10.1021/ja054872i and 10.1021/ja00311a085 and

10.1021/ja00122a011

Retrosynthesis ID: 9999983

2.5 Path 5

Score: 31.25




Figure 5: Outline of path 5

2.5.1 Suzuki coupling of alkyl-9-BBNs with vinyl iodides

Substrates:

1. 1-bromo-3-(2-iodoethenyl)benzene - available at Sigma-Aldrich

2. 9-(3,3-pentamethyleneallyl)-9-borabicyclo3.3.1nonane

Products:

1. Brc1cccc(/C=C/CC=C2CCCC2)c1

 ${\bf Typical\ conditions:}\ {\bf Pd\ catalyst.base.solvent}$

Protections: none

Reference: 10.1021/jo015995y and 10.1016/j.tetlet.2010.11.139 And 10.1021/ol0600741 and 10.1055/s-2002-32602 and 10.1002/anie.200501760

Retrosynthesis ID: 25168

2.5.2 Br/Li exchange

1. n-BuLi - available at Sigma-Aldrich

 $2. \ \, Brc1cccc(/C=C/CC=C2CCCC2)c1$

Products:

1. [Li]c1cccc(/C=C/CC=C2CCCC2)c1


Typical conditions: nBuLi.or.tBuLi.THF.-78C

Protections: none

Reference: 10.1002/ejoc.201101490 and 10.1016/j.tet.2012.03.058 and 10.1016/j.tetlet.2015.01.032 and 10.1021/ja0541175 and 10.1016/j.tetlet.2016.06.123

Retrosynthesis ID: 30672

2.5.3 Addition of electrophiles to lithiated arenes/heteroarenes



Substrates:

1. [Li]c1cccc(/C=C/CC=C2CCCC2)c1

Products:

1. Oc1cccc(/C=C/CC=C2CCCC2)c1

Typical conditions: B(OMe)3 then H2O2.THF

Protections: none

Reference: 10.1039/C7CC09187G (SI) and 10.1002/ejoc.201701142 and

10.1021/acscatal.6b03380 (SI,p.10) and 10.1002/chem.201702143 (SI)