

Paths of analysis*

L8_DIA

Synthia

October 11, 2022

1 Analysis parameters

Analysis type: Automatic Retrosynthesis

Rules: none selected

Filters: Tunnels, FGI, FGI with protections

Max. paths returned: 50

Max. iterations: 2000

Commercial:

1. Max. molecular weight - 1000 g/mol
2. Max. price - 1500 \$/g

Published:

1. Max. molecular weight - 1000 g/mol
2. Popularity - 5

My Stockroom:

1. Max. molecular weight - 1000 g/mol

Reaction scoring formula: $\text{TUNNEL_COEF} * \text{FGI_COEF} * \text{STEP} * 20 + 1000000 * (\text{CONFLICT} + \text{NON_SELECTIVITY} + \text{FILTERS} + \text{PROTECT})$

Chemical scoring formula: $\text{SMALLER}^3, \text{SMALLER}^{1.5}$

Min. search width: 400

Max. reactions per product: 60

Strategies: none selected

*The results stated herein were generated using the proprietary platform owned and maintained by Grzybowski Scientific Inventions, Inc., a subsidiary of Merck KGaA, Darmstadt Germany. The results are provided on an as is basis, and shall be used solely in connection with the rights afforded in the license agreement and for no other purpose.

FGI Coeff: 0

Tunnels Coeff: 0

JSON Parameters: {}

2 Paths

4 paths found. *Paths are sorted by score. Reactions are sorted in appearance order for each path.*

2.1 Path 1

Score: 100.08





Figure 1: Outline of path 1

2.1.1 Steglich Esterification



Substrates:

1. crotonalkohol

2. 5-phenyl-penta-2,4-dienoic acid

Products:

1. CC=CCOC(=O)C=CC=Cc1ccccc1


Typical conditions: alcohol.DCC.DMAP.DCM or thiol.DCC.DMAP.DCM

Protections: none

Reference: [10.1002/anie.197805221](#)

Retrosynthesis ID: 10171

2.1.2 Wohl-Ziegler Bromination



Substrates:

1. N-Bromosuccinimide - [available at Sigma-Aldrich](#)
2. CC=CCOC(=O)C=CC=Cc1ccccc1

Products:

1. O=C(C=CC=Cc1ccccc1)OCC=CCBr


Typical conditions: NBS.AIBN or (BzO)₂ or heat

Protections: none

Reference: [10.1016/j.steroids.2018.10.005](#) (Scheme 1) and [10.1016/j.bmc.2010.06.075](#) (Scheme 2) and [10.1021/acs.orglett.9b03865](#) (p. SI 6)

Retrosynthesis ID: 245554

2.1.3 Diels-Alder



Substrates:

1. O=C(C=CC=CC1CCCCC1)OCC=CCBr

Products:

1. O=C1OCC2C1C=CC(c1cccc1)C2CBr


Typical conditions: Lewis acid or chiral Lewis acid. Solvent.

Protections: none

Reference: DOI: [10.1002/1521-3773\(20020517\)41:10<1668::AID-ANIE1668>3.0.CO;2-Z](https://doi.org/10.1002/1521-3773(20020517)41:10<1668::AID-ANIE1668>3.0.CO;2-Z) AND [10.1021/ja062508t](https://doi.org/10.1021/ja062508t)

Retrosynthesis ID: 18116

2.1.4 Wittig-Schlosser olefination



Substrates:

1. O=C1OCC2C1C=CC(c1ccccc1)C2CBr
2. Ethanal - *available at Sigma-Aldrich*

Products:


1. C/C=C/C1C(c2ccccc2)C=CC2C(=O)OCC21

Typical conditions: 1.PPh₃ or trialkylphosphite.2.base.aldehyde.3.base

Protections: none

Reference: [10.1021/ol049701h](#) and [10.1021/ja00535a063](#) and Kurti and Czako; Strategic Applications of Named Reactions in Organic Synthesis. 1st edn., 488-489.

Retrosynthesis ID: 9546

2.1.5 Claisen Condensation**Substrates:**

1. C/C=C/C1C(c2ccccc2)C=CC2C(=O)OCC21
2. Methyl acetate - *available at Sigma-Aldrich*

Products:

1. C/C=C/C1C(c2ccccc2)C=CC2(C(C)=O)C(=O)OCC12

Typical conditions: Base.Solvent

Protections: none

Reference: [10.1021/cr020703u](#) and [10.1021/cr60088a002](#)

Retrosynthesis ID: 5015

2.2 Path 2

Score: 105.08





Figure 2: Outline of path 2

2.2.1 Takai olefination



Substrates:

1. butyryloxy-acetaldehyde
2. Bromoform - *available at Sigma-Aldrich*

Products:

1. CCCC(=O)OCC=CCBr


Typical conditions: CrCl₂.THF

Protections: none

Reference: [10.1021/ja00283a046](#) and [10.1021/ja00237a081](#)

Retrosynthesis ID: 11463

2.2.2 Condensation of esters with aldehydes/ketones



Substrates:

1. cinnamaldehyd - *available at Sigma-Aldrich*
2. CCCC(=O)OCC=CBr

Products:

1. CCC(=CC=Cc1ccccc1)C(=O)OCC=CBr


Typical conditions: LDA.THF

Protections: none

Reference: [10.1021/op040006z](#) AND [10.1016/j.bmcl.2005.10.104](#) AND

Retrosynthesis ID: 14983

2.2.3 Allylic Oxidation of Alkenes



Substrates:

1. CCC(=CC=Cc1ccccc1)C(=O)OCC=CBr

Products:

1. CC(=O)C(=CC=Cc1ccccc1)C(=O)OCC=CBr


Typical conditions: tBuOOH.Pd(OH)2/C or PhI(OAc)2 or SeO2

Protections: none

Reference: [10.1021/ja0340735](#) and [10.1021/ol100603q](#) and [10.1016/j.tetlet.2016.05.063](#) (Scheme 2)

Retrosynthesis ID: 2583

2.2.4 Pd-catalyzed formylation of vinyl halides



Substrates:

1. tert-Butyl isocyanide - *available at Sigma-Aldrich*
2. CC(=O)C(=CC=Cc1ccccc1)C(=O)OCC=CBr

Products:

1. CC(=O)C(=CC=Cc1ccccc1)C(=O)OCC=CC=O


Typical conditions: Pd(OAc)₂.JohnPhos.Na₂CO₃.H₂O.Et₃SiH.DMF.65°C

Protections: none

Reference: DOI: [10.1021/ol5014262](https://doi.org/10.1021/ol5014262)

Retrosynthesis ID: 3104

2.2.5 Diels-Alder



Substrates:

1. CC(=O)C(=CC=Cc1ccccc1)C(=O)OCC=CC=O

Products:

1. CC(=O)C12C=CC(c3ccccc3)C(C=O)C1COC2=O


Typical conditions: Lewis acid or chiral Lewis acid. Solvent.

Protections: none

Reference: DOI: [10.1002/1521-3773\(20020517\)41:10<1668::AID-ANIE1668>3.0.CO;2-Z](https://doi.org/10.1002/1521-3773(20020517)41:10<1668::AID-ANIE1668>3.0.CO;2-Z) AND [10.1021/ja062508t](https://doi.org/10.1021/ja062508t)

Retrosynthesis ID: 18116

2.2.6 Takai olefination



Substrates:

- 1,1-Diiodoethane - *available at Sigma-Aldrich*
- CC(=O)C12C=CC(c3ccccc3)C(C=O)C1COC2=O

Products:

- C/C=C/C1C(c2ccccc2)C=CC2(C(C)=O)C(=O)OCC12

Typical conditions: CrCl₂.THF.DMF

Protections: none

Reference: [10.1021/ja00283a046](#) and [10.1021/ja00237a081](#)

Retrosynthesis ID: 10942

2.3 Path 3

Score: 120.37





Figure 3: Outline of path 3

2.3.1 Brown Hydroboration of Alkenes

**Substrates:**

- Allyl butyrate - *available at Sigma-Aldrich*

Products:

- butyric acid 3-hydroxy-propyl ester


Typical conditions: B₂H₆.H₂O₂.THF.NaOH

Protections: none

Reference: [10.1002/9780470638859.conrr118](https://doi.org/10.1002/9780470638859.conrr118)

Retrosynthesis ID: 4772

2.3.2 Oxidation of primary alcohols with DMP



Substrates:

1. butyric acid 3-hydroxy-propyl ester

Products:

1. CCCC(=O)OCCC=O


Typical conditions: DMP.DCM.0-25 C

Protections: none

Reference: [10.1016/j.bmc.2020.115469](https://doi.org/10.1016/j.bmc.2020.115469) p. 3, 9 and [10.1021/acs.jmedchem.8b01878](https://doi.org/10.1021/acs.jmedchem.8b01878) SI p. S43

Retrosynthesis ID: 50426

2.3.3 Tebbe Olefination



Substrates:

1. CCCC(=O)OCCC=O

Products:

1. 3-butenyl butyrate


Typical conditions: Cp2TiCl2.AlMe3.toluene

Protections: none

Reference: [10.1016/j.tet.2007.03.015](https://doi.org/10.1016/j.tet.2007.03.015) and [10.1002/9780470638859.conrr617](https://doi.org/10.1002/9780470638859.conrr617)

Retrosynthesis ID: 11714

2.3.4 Condensation of esters with aldehydes/ketones



Substrates:

1. 3-butenyl butyrate
2. cinnamaldehyde - *available at Sigma-Aldrich*

Products:

1. C=CCCOC(=O)C(=CC=Cc1ccccc1)CC


Typical conditions: LDA.THF

Protections: none

Reference: [10.1021/op040006z](#) AND [10.1016/j.bmcl.2005.10.104](#) AND

Retrosynthesis ID: 14983

2.3.5 Tsuji-Wacker Oxidation of alkenes



Substrates:

1. C=CCCOC(=O)C(=CC=Cc1ccccc1)CC

Products:

1. CCC(=CC=Cc1ccccc1)C(=O)OCCCC=O


Typical conditions: PdCl₂(PhCn)₂.CuCl₂.AgNO₂.O₂.tBuOH.MeNO₂.rt

Protections: none

Reference: [10.1021/jacs.6b08788](#) and [10.1021/ja411749k](#) and [10.1002/anie.201306756](#) and [10.1016/S0040-4039\(03\)01709-X](#) and [10.1021/acs.orglett.6b01165](#)

Retrosynthesis ID: 28273

2.3.6 Enol esters and ethers synthesis



Substrates:

1. CCC(=CC=Cc1ccccc1)C(=O)OCCCC=O
2. TMSCl - [available at Sigma-Aldrich](#)

Products:

1. CCC(=CC=Cc1ccccc1)C(=O)OCCC=CO[Si](C)(C)C


Typical conditions: 1.LDA.2.Electrophile

Protections: none

Reference: US2467095A AND WO2014169833a1 AND [10.1016/j.steroids.2011.03.014](#) AND [10.1021/ol200875m](#) (SI) AND [10.1021/ja00531a034](#)

Retrosynthesis ID: 7797

2.3.7 Dehydrogenation of silyl enol ethers



Substrates:

1. CCC(=CC=Cc1ccccc1)C(=O)OCCC=CO[Si](C)(C)C

Products:

1. CCC(=CC=Cc1ccccc1)C(=O)OCC=CC=O


Typical conditions: Pd(OAc)₂.Cu(OAc)₂.O₂.MeCN

Protections: none

Reference: [10.1271/bbb.60.405](#) and [10.1039/C3CC46778C](#) and US2015284405 p.40 and [10.1016/S0040-4039\(01\)81518-5](#) and US2010204477 p. 15-16 and [10.1016/0040-4039\(95\)00694-8](#) and [10.1021/jo00089a034](#) and [10.1016/S0040-4020\(01\)90587-3](#) and [10.1080/00397919008052802](#) and [10.1021/ja00218a060](#)

Retrosynthesis ID: 9999877

2.3.8 Reformatsky Reaction



Substrates:

1. Methyl 2-bromopropionate - *available at Sigma-Aldrich*
2. CCC(=CC=Cc1ccccc1)C(=O)OCC=CC=O

Products:

1. CCC(=CC=Cc1ccccc1)C(=O)OCC=CC(O)C(C)C(=O)OC


Typical conditions: Me₂Zn.B(OMe)₃.toluene.Et₂O

Protections: none

Reference: [10.1021/jo200774e](#) p. 6373 and [10.1021/jo00163a019](#) p. 2522, 2525

Retrosynthesis ID: 11164

2.3.9 Allylic Oxidation of Alkenes



Substrates:

1. CCC(=CC=Cc1ccccc1)C(=O)OCC=CC(O)C(C)C(=O)OC

Products:

1. COC(=O)C(C)C(O)C=CCOC(=O)C(=CC=Cc1ccccc1)C(C)=O


Typical conditions: tBuOOH.Pd(OH)₂/C or PhI(OAc)₂ or SeO₂

Protections: none

Reference: [10.1021/ja0340735](#) and [10.1021/ol100603q](#) and [10.1016/j.tetlet.2016.05.063](#) (Scheme 2)

Retrosynthesis ID: 2583

2.3.10 Diels-Alder



Substrates:

1. COC(=O)C(C)C(O)C=CCOC(=O)C(=CC=Cc1ccccc1)C(C)=O

Products:

1. COC(=O)C(C)C(O)C1C(c2ccccc2)C=CC2(C(C)=O)C(=O)OCC12


Typical conditions: Lewis acid or chiral Lewis acid. Solvent.

Protections: none

Reference: DOI: [10.1002/1521-3773\(20020517\)41:10<1668::AID-ANIE1668>3.0.CO;2-Z](https://doi.org/10.1002/1521-3773(20020517)41:10<1668::AID-ANIE1668>3.0.CO;2-Z) AND [10.1021/ja062508t](https://doi.org/10.1021/ja062508t)

Retrosynthesis ID: 18116

2.3.11 Tandem Krapcho decarboxylation and elimination



Substrates:

1. glacial - *available at Sigma-Aldrich*
2. COC(=O)C(C)C(O)C1C(c2ccccc2)C=CC2(C(C)=O)C(=O)OCC12

Products:

1. C/C=C/C1C(c2ccccc2)C=CC2(C(C)=O)C(=O)OCC12

Typical conditions: 1. Ac₂O.py 2. DMSO.H₂O.NaCl.170C

Protections: none

Reference: DOI: [10.1021/jo00263a005](https://doi.org/10.1021/jo00263a005) and [10.1021/jo00386a011](https://doi.org/10.1021/jo00386a011) and [10.1021/ol006085q](https://doi.org/10.1021/ol006085q)

Retrosynthesis ID: 9605

2.4 Path 4

Score: 129.38





Figure 4: Outline of path 4

2.4.1 Enol esters and ethers synthesis



Substrates:

1. TMSCl - *available at Sigma-Aldrich*
2. hex-4t-enal

Products:

1. C/C=C/CC=CO[Si](C)(C)C


Typical conditions: 1.LDA.2.Electrophile

Protections: none

Reference: US2467095A AND WO2014169833a1 AND
[10.1016/j.steroids.2011.03.014](#) AND [10.1021/ol200875m](#) (SI) AND
[10.1021/ja00531a034](#)

Retrosynthesis ID: 7797

2.4.2 Dehydrogenation of silyl enol ethers



Substrates:

1. C/C=C/CC=CO[Si](C)(C)C

Products:

1. 2,4-hexadienal


Typical conditions: Pd(OAc)₂.Cu(OAc)₂.O₂.MeCN

Protections: none

Reference: [10.1271/bbb.60.405](#) and [10.1039/C3CC46778C](#) and US2015284405 p.40 and [10.1016/S0040-4039\(01\)81518-5](#) and US2010204477 p. 15-16 and [10.1016/0040-4039\(95\)00694-8](#) and [10.1021/jo00089a034](#) and [10.1016/S0040-4020\(01\)90587-3](#) and [10.1080/00397919008052802](#) and [10.1021/ja00218a060](#)

Retrosynthesis ID: 9999877

2.4.3 Synthesis of alkyl chlorides from aldehydes



Substrates:

1. 2,4-hexadienal

Products:

1. C/C=C/C=CCCCl


Typical conditions: InO3.chloroform.SiMe2Cl

Protections: none

Reference: DOI: [10.1021/ja0283246](https://doi.org/10.1021/ja0283246)

Retrosynthesis ID: 11621

2.4.4 Aldol Condensation



Substrates:

1. 4-Hydroxy-2-butanone - *available at Sigma-Aldrich*
2. cinnamaldehyde - *available at Sigma-Aldrich*

Products:

1. CC(=O)C(=CC=Cc1ccccc1)CO


Typical conditions: NaOEt.base

Protections: none

Reference: [10.1080/00397911.2016.1206938](#)

Retrosynthesis ID: 10049

2.4.5 Diels-Alder



Substrates:

1. C/C=C/C=C\CCl
2. CC(=O)C(=CC=Cc1ccccc1)CO

Products:

1. C/C=C/C1C(c2ccccc2)C=CC(CO)(C(C)=O)C1CCl


Typical conditions: Lewis acid or chiral Lewis acid. Solvent.

Protections: none


Reference: DOI: [10.1002/1521-3773\(20020517\)41:10<1668::AID-ANIE1668>3.0.CO;2-Z](#) AND [10.1021/ja062508t](#)

Retrosynthesis ID: 18116


2.4.6 Jones Oxidation



Substrates:



Products:




Typical conditions: cromate.sulfate.H2O.acetone

Protections: none

Reference: [10.1002/9780470638859.conrr349](#) and [10.1021/jm00270a004](#)

Retrosynthesis ID: 11160

2.4.7 Synthesis of esters from alkyl chlorides and carboxylic acids or thioacids



Substrates:

1. C/C=C/C1C(c2ccccc2)C=CC(C(C)=O)(C(=O)O)C1CCl

Products:

1. C/C=C/C1C(c2ccccc2)C=CC2(C(C)=O)C(=O)OCC12

Typical conditions: K₂CO₃.DMF

Protections: none

Reference: [10.1016/j.bmcl.2005.08.026](#) AND [10.1021/ol034655r](#) (SI) AND [10.1039/C3RA41967C](#) AND [10.1016/j.bmcl.2012.03.093](#)

Retrosynthesis ID: 14684