Paths of analysis*

Synthia

October 10, 2022

1 Analysis parameters

Analysis type: Automatic Retrosynthesis

Rules: none selected

Filters: Exclude Diastereoselecitve reactions, Tunnels, FGI, FGI with protec-

tions

Max. paths returned: 50

Max. iterations: 2000

Commercial:

1. Max. molecular weight - 1000 g/mol

2. Max. price - 1500 \$/g

Published:

- 1. Max. molecular weight 1000 g/mol
- 2. Popularity 5

My Stockroom:

1. Max. molecular weight - 1000 g/mol

Reaction scoring formula: TUNNEL_COEF*FGI_COEF*STEP*20+1000 000*(CONFLICT+NON SELECTIVITY+FILTERS+PROTECT)

Chemical scoring formula: SMALLER^ 3,SMALLER^ 1.5

Min. search width: 400

Max. reactions per product: 60

^{*}The results stated herein were generated using the proprietary platform owned and maintained by Grzybowski Scientific Inventions, Inc., a subsidiary of Merck KGaA, Darmstadt Germany. The results are provided on an as is basis, and shall be used solely in connection with the rights afforded in the license agreement and for no other purpose.

Strategies: none selected

FGI Coeff: 0

Tunnels Coeff: 0

JSON Parameters: {}

2 Paths

5 paths found. Paths are sorted by score. Reactions are sorted in appearance order for each path.


2.1 Path 1

Score: 125.08



Figure 1: Outline of path 1

2.1.1 Oxidation of primary alcohols with DMP



Substrates:

1. 2-phenylbut-3-en-1-ol

Products:

1. 2-phenyl-but-3-enal

Typical conditions: DMP.DCM.0-25 $\rm C$


Protections: none

Reference: 10.1016/j.bmc.2020.115469 p. 3, 9 and

10.1021/acs.jmedchem.8b01878 SI p. S43

Retrosynthesis ID: 50426

2.1.2 Iodoolefination of aldehydes



Substrates:

1. Diiodomethane - available at Sigma-Aldrich

2. 2-phenyl-but-3-enal

Products:

1. $C=CC(/C=C\setminus I)c1ccccc1$


Typical conditions: 1.PPh3.2.NaN(TMS)2.HMPA.THF

Protections: none

Reference: 10.1021/ja00171a035 and 10.1039/C0OB00977F and WO2009033499

(p.25)

2.1.3 Rh(I) catalyzed hydroboration



Substrates:

1. $C=CC(/C=C\setminus I)c1ccccc1$

Products:

1. $OCCC(/C=C\setminus I)c1ccccc1$


 $\textbf{Typical conditions:} \ \ \textbf{Wilkinson's catalyst.} cateholborane. THF. MeOH. NaOH. H2O2$

Protections: none

Reference: DOI: 10.1021/ja00043a009

Retrosynthesis ID: 9910000

2.1.4 Stille Reaction of Vinyl Iodides with Alkenyl Stannanes



Substrates:

1. $OCCC(/C=C\backslash I)c1ccccc1$

2. 3-tributylstannyl-2(5h)-furanone

Products:

1. $O=C1OCC=C1/C=C\setminus C(CCO)c1ccccc1$


Typical conditions: Pd (cat). Ligand. CuCl. THF. AsPh3. Heating

Protections: none

Reference: DOI: 10.1021/jo047732k or 10.1021/ja062524q (SI, page S35) or 10.1021/ja1009579 (SI, page S-15) or 10.1025/s-0034-1378360 or 10.1021/ja028726d or 10.1016/j.tetlet.2015.10.087

Retrosynthesis ID: 9991429

2.1.5 Oxidation of primary alcohols with DMP



Substrates:

1. $O=C1OCC=C1/C=C\setminus C(CCO)c1ccccc1$

Products:

1. $O=CCC(/C=C\setminus C1=CCOC1=O)c1ccccc1$

Typical conditions: DMP.DCM.0-25 C

Protections: none

Reference: 10.1016/j.bmc.2020.115469 p. 3, 9 and 10.1021/acs.jmedchem.8b01878 SI p. S43

2.1.6 Michael addition

Substrates:

1. $O=CCC(/C=C\setminus C1=CCOC1=O)c1ccccc1$

Products:

1. O=CC1C(c2cccc2)C=CC2C(=O)OCC21


 $\textbf{Typical conditions:} \ \, \textbf{EtONa or other base}$

Protections: none

Reference: 10.1016/j.tetlet.2011.02.073 AND 10.1016/j.molstruc.2010.12.005 AND 10.1016/S0040-4039(97)00695-3 AND 10.1021/ol016401g AND 10.1002/ejoc.200500330

Retrosynthesis ID: 15774

2.1.7 Wittig-Schlosser olefination



Substrates:

1. Bromoethane - available at Sigma-Aldrich

 $2. \ O=CC1C(c2cccc2)C=CC2C(=O)OCC21$

Products:

1. C/C=C/C1C(c2cccc2)C=CC2C(=O)OCC21

Typical conditions: 1.PPh3 or trialkylphosphite.2.base.aldehyde.3.base

Protections: none

Reference: 10.1021/ol049701h and 10.1021/ja00535a063 and Kurti and Czako; Strategic Applications of Named Reactions in Organic Synthesis. 1st edn., 488-489.

Retrosynthesis ID: 9546

2.1.8 Claisen Condensation

Substrates:

1. C/C=C/C1C(c2cccc2)C=CC2C(=O)OCC21

2. Methyl acetate - available at Sigma-Aldrich

Products:

 $1. \ C/C = C/C1C(c2cccc2)C = CC2(C(C) = O)C(=O)OCC12$

Typical conditions: Base.Solvent

Protections: none

Reference: 10.1021/cr020703u and 10.1021/cr60088a002

Retrosynthesis ID: 5015

2.2 Path 2

Score: 146.56





Figure 2: Outline of path 2

2.2.1 Coupling of alkynes and alcohols



${\bf Substrates:}$

1. 2-Buten-1-ol - available at Sigma-Aldrich

2. 1-Phenyl-1-propyne - available at Sigma-Aldrich

Products:

1. C13H16O

Typical conditions: H2Ru(CO)(PPh3)3.2,4,6-(iPr)3PhSO3H.SL-J009-

1.TBAI.IPA.THF.95C

Protections: none

Reference: DOI: 10.1021/jacs.5b00747

2.2.2 Appel Reaction

Substrates:

1. C13H16O

Products:

 $1. \ C{=}CC(c1ccccc1)C(Br)/C{=}C/C$


Typical conditions: PPh3.CBr4

 ${\bf Protections:}\ {\rm none}$

Reference: 10.1016/j.jfluchem.2015.03.009 and 10.1016/j.tet.2005.12.006 and 10.1021/jm00161a029 and 10.1055/s-1995-5215

Retrosynthesis ID: 9990042

2.2.3 Suzuki coupling of vinyl bromides with alkenyl boronic acids



Substrates:

1. Vinylboronic acid - available at Sigma-Aldrich

2. 3-bromo-2,5-dihydrofuran-2-one - available at Sigma-Aldrich

Products:

1. 3-vinyl-2(5h)-furanone

Typical conditions: Pd catalyst.base.solvent

Protections: none

Reference: 10.1021/cr00039a007 and $10.1007/3418_2012_32$ and 10.1021/cr0505268 and 10.1016/j.jfluchem.2016.01.018 and 10.1039/C3CS60197H

Retrosynthesis ID: 24926

2.2.4 Conjugated addition of organocuprate-acylation of enones and enoate esters

Substrates:

- 1. C=CC(c1cccc1)C(Br)/C=C/C
- 2. 3-vinyl-2(5h)-furanone
- 3. Acetyl chloride available at Sigma-Aldrich

Products:

 $1. \ C = CC(c1ccccc1)C(/C = C/C)C1COC(=O)C1(C = C)C(C) = O$

Typical conditions: 1.RCuLi.2.AcCl.HMPA

Protections: none

Reference: 10.3987/COM-99-S143 AND 10.1021/ja00148a023 AND

10.1016/S0040-4039(01)80891-1

${\bf 2.2.5} \quad {\bf Ring\text{-}Closing\ Metathesis}$

Substrates:

 $1. \ C{=}CC(c1ccccc1)C(/C{=}C/C)C1COC({=}O)C1(C{=}C)C(C){=}O\\$

Products:

 $1. \ C/C = C/C1C(c2cccc2)C = CC2(C(C) = O)C(=O)OCC12$

 $\textbf{Typical conditions:} \ \ \text{catalyst e.g.} \ \ \text{Hoveyda-Grubbs} \ \ . \ \ \text{solvent e.g.} \ \ \text{CH2Cl2}$

Protections: none

Reference: DOI: 10.1002/anie.200800693 and 10.1021/acs.orglett.8b04003 and 10.1021/jo0264729 and 10.1021/ja072334v and 10.1002/ejoc.201001102

Retrosynthesis ID: 31014187

2.3 Path 3

Score: 161.11




Figure 3: Outline of path 3

2.3.1 Corey-Seebach

Substrates:

- 1. (e)-1-(1,3-dithian-2-yl)but-2-ene
- 2. Formalin available at Sigma-Aldrich

Products:

1. C/C=C/CC(=O)CO

Typical conditions: BuLi.THF.-30C.HgO.H2O.THF

Protections: none

Reference: 10.1055/s-1977-24412

Retrosynthesis ID: 11019

2.3.2 Robinson annulation

Substrates:

$$1. \ \mathrm{C/C}{=}\mathrm{C/CC}(=\mathrm{O})\mathrm{CO}$$

2. 4-phenylbutenone

Products:

1. C/C=C/C1C(CO)=CC(=O)CC1c1ccccc1

Typical conditions: KOH.MeOH

Protections: none

Reference: DOI: 10.1021/ja9602509 (SI) AND DOI: 10.1021/ja00735a059 AND

10.1021/ol1011955

Retrosynthesis ID: 7595

2.3.3 Lactonization of allylic alcohols

Substrates:

1. CORM-2 - available at Sigma-Aldrich

2. C/C=C/C1C(CO)=CC(=O)CC1c1ccccc1

Products:

 $1. \ C/C=C/C1C(c2cccc2)CC(=O)C2C(=O)OCC21$

 $\textbf{Typical conditions:} \ \ O2.PdCl2.CuCl2.HCl.THF.rt$

Protections: none

Reference: DOI: 10.1016/S0040-4039(01)80907-2 and 10.1021/jo9702709

2.3.4 Synthesis of enol phosphonates

Substrates:

- $1. \ \mathrm{C/C}{=}\mathrm{C/C1C}(\mathrm{c2cccc2})\mathrm{CC}({=}\mathrm{O})\mathrm{C2C}({=}\mathrm{O})\mathrm{OCC21}$
- 2. Diethyl chlorophosphate available at Sigma-Aldrich

Products:

 $1. \ \ C/C=C/C1C(c2cccc2)C=C(OP(=O)(OCC)OCC)C2C(=O)OCC21$

Typical conditions: 1. LDA.2.ClP(=O)(OR)2

Protections: none

Reference: DOI: 10.1016/j.tetasy.2012.09.012

Retrosynthesis ID: 7573

2.3.5 Claisen Condensation

Substrates:

 $1. \ C/C=C/C1C(c2cccc2)C=C(OP(=O)(OCC)OCC)C2C(=O)OCC21$

2. Methyl acetate - available at Sigma-Aldrich

Products:

 $1. \ C/C = C/C1C(c2cccc2)C = C(OP(=O)(OCC)OCC)C2(C(C)=O)C(=O)OCC12$


Typical conditions: Base.Solvent

Protections: none

Reference: 10.1021/cr020703u and 10.1021/cr60088a002

Retrosynthesis ID: 5015

2.3.6 Reduction of enol phosphonates



Substrates:

 $1. \ C/C = C/C1C(c2cccc2)C = C(OP(=O)(OCC)OCC)C2(C(C)=O)C(=O)OCC12$

Products:

 $1. \ C/C = C/C1C(c2cccc2)C = CC2(C(C) = O)C(=O)OCC12$

Typical conditions: Et3Al.Pd(PPh3)4

Protections: none

Reference: 10.1021/jo00387a038 AND 10.1021/jo00292a049 AND

10.1039/C1CS15100B

2.4 Path 4

Score: 164.14





Figure 4: Outline of path 4

2.4.1 Coupling of alkynes and alcohols



Substrates:

1. 2-Buten-1-ol - available at Sigma-Aldrich

2. 1-Phenyl-1-propyne - available at Sigma-Aldrich

Products:

1. C13H16O

Typical conditions: H2Ru(CO)(PPh3)3.2,4,6-(iPr)3PhSO3H.SL-J009-

 $1.\mathrm{TBAI.IPA.THF.95C}$

Protections: none

Reference: DOI: 10.1021/jacs.5b00747

2.4.2 Sulfonation of secondary alcohols

Substrates:

- 1. C13H16O
- 2. Mesyl chloride available at Sigma-Aldrich

Products:

1. C=CC(c1cccc1)C(/C=C/C)OS(C)(=O)=O

 $\textbf{Typical conditions:} \ Et 3N.DMAP.DCM$

 ${\bf Protections:}\ {\rm none}$

Reference: 10.1021/jo048289g and 10.1021/ja9617808 and

10.1016/j.steroids.2005.10.004

Retrosynthesis ID: 24386

2.4.3 Substitution of secondary mesyl group with bromide



Substrates:

1. C=CC(c1cccc1)C(/C=C/C)OS(C)(=O)=O

Products:

1. C=CC(c1cccc1)C(Br)/C=C/C


Typical conditions: LiBr.DMF

Protections: none

Reference: 10.1021/jo00068a037 and 10.1016/S0040-4020(03)00140-6

Retrosynthesis ID: 29713

2.4.4 Suzuki coupling of vinyl bromides with alkenyl boronic acids



Substrates:

- 1. Vinylboronic acid available at Sigma-Aldrich
- 2. 3-bromo-2,5-dihydrofuran-2-one available at Sigma-Aldrich

Products:

1. 3-vinyl-2(5h)-furanone

Typical conditions: Pd catalyst.base.solvent

Protections: none

Reference: 10.1021/cr00039a007 and $10.1007/3418_2012_32$ and 10.1021/cr0505268 and 10.1016/j.jfluchem.2016.01.018 and 10.1039/C3CS60197H

2.4.5 Conjugated addition of organocuprate-acylation of enones and enoate esters

Substrates:

- 1. C=CC(c1cccc1)C(Br)/C=C/C
- 2. 3-vinyl-2(5h)-furanone
- 3. Acetyl chloride available at Sigma-Aldrich

Products:

1. C=CC(c1ccccc1)C(/C=C/C)C1COC(=O)C1(C=C)C(C)=O

Typical conditions: 1.RCuLi.2.AcCl.HMPA

Protections: none

Reference: 10.3987/COM-99-S143 AND 10.1021/ja00148a023 AND

10.1016/S0040-4039(01)80891-1

Retrosynthesis ID: 20523

2.4.6 Ring-Closing Metathesis

Substrates:

 $1. \ C = CC(c1ccccc1)C(/C = C/C)C1COC(=O)C1(C = C)C(C) = O$

Products:

1. C/C=C/C1C(c2cccc2)C=CC2(C(C)=O)C(=O)OCC12

Typical conditions: catalyst e.g. Hoveyda-Grubbs . solvent e.g. CH2Cl2

Protections: none

 $\textbf{Reference:} \ \ DOI: \ \textit{10.1002/anie.200800693} \ \ \text{and} \ \ \textit{10.1021/acs.orglett.8b04003} \ \ \text{and}$

10.1021/jo0264729 and 10.1021/ja072334v and 10.1002/ejoc.201001102

Retrosynthesis ID: 31014187

2.5 Path 5

Score: 164.14





Figure 5: Outline of path 5

2.5.1 Suzuki coupling of vinyl bromides with alkenyl boronic acids



Substrates:

- 1. Vinylboronic acid available at Sigma-Aldrich
- 2. 3-bromo-2,5-dihydrofuran-2-one available at Sigma-Aldrich

Products:

1. 3-vinyl-2(5h)-furanone

 ${\bf Typical\ conditions:}\ {\bf Pd\ catalyst.base.solvent}$

Protections: none

Reference: 10.1021/cr00039a007 and 10.1007/3418_2012_32 and 10.1021/cr0505268 and 10.1016/j.jfluchem.2016.01.018 and 10.1039/C3CS60197H

Retrosynthesis ID: 24926

2.5.2 Michael addition

Substrates:

- 1. 3-vinyl-2(5h)-furanone
- 2. 3-phenyl-4-pentenal

Products:

 $1. \ C{=}CC1C({=}O)OCC1C(C{=}O)C(C{=}C)c1ccccc1$

 $\textbf{Typical conditions:} \ \, \textbf{EtONa or other base}$

Protections: none

Reference: 10.1016/j.tetlet.2011.02.073 AND 10.1016/j.molstruc.2010.12.005 AND 10.1016/S0040-4039(97)00695-3 AND 10.1021/ol016401g AND 10.1002/ejoc.200500330

2.5.3 Takai olefination

Substrates:

1. 1,1-Diiodoethane - available at Sigma-Aldrich

 $2. \ C{=}CC1C({=}O)OCC1C(C{=}O)C(C{=}C)c1ccccc1$

Products:

1. C=CC1C(=O)OCC1C(/C=C/C)C(C=C)c1ccccc1


Typical conditions: CrCl2.THF.DMF

Protections: none

Reference: 10.1021/ja00283a046 and 10.1021/ja00237a081

Retrosynthesis ID: 10942

2.5.4 Claisen Condensation



Substrates:

1. Methyl acetate - available at Sigma-Aldrich

 $2. \ C{=}CC1C({=}O)OCC1C(/C{=}C/C)C(C{=}C)c1ccccc1\\$

Products:

 $1. \ C{=}CC(c1ccccc1)C(/C{=}C/C)C1COC({=}O)C1(C{=}C)C(C){=}O\\$

 ${\bf Typical\ conditions:}\ {\bf Base. Solvent}$

Protections: none

Reference: 10.1021/cr020703u and 10.1021/cr60088a002

Retrosynthesis ID: 5015

2.5.5 Ring-Closing Metathesis

Substrates:

 $1. \ C = CC(c1ccccc1)C(/C = C/C)C1COC(=O)C1(C = C)C(C) = O$

Products:

 $1. \ C/C = C/C1C(c2cccc2)C = CC2(C(C) = O)C(=O)OCC12$

 $\textbf{Typical conditions:} \ \, \text{catalyst e.g. Hoveyda-Grubbs} \,\, . \,\, \text{solvent e.g. CH2Cl2}$

Protections: none

Reference: DOI: 10.1002/anie.200800693 and 10.1021/acs.orglett.8b04003 and

 $10.1021/jo0264729 \ \ {\rm and} \quad 10.1021/ja072334v \ \ {\rm and} \quad 10.1002/ejoc.201001102$