Paths of analysis*

Synthia

October 10, 2022

1 Analysis parameters

Analysis type: Automatic Retrosynthesis

Rules: none selected

Filters: Exclude Diastereoselective reactions, Tunnels, FGI, FGI with protec-

tions

Max. paths returned: 50

Max. iterations: 2000

Commercial:

1. Max. molecular weight - 1000 g/mol

2. Max. price - 1500 \$/g

Published:

- 1. Max. molecular weight 1000 g/mol
- 2. Popularity 5

My Stockroom:

1. Max. molecular weight - 1000 g/mol

Reaction scoring formula: TUNNEL_COEF*FGI_COEF*STEP*20+1000 000*(CONFLICT+NON SELECTIVITY+FILTERS+PROTECT)

Chemical scoring formula: SMALLER^ 3,SMALLER^ 1.5

Min. search width: 400

Max. reactions per product: 60

^{*}The results stated herein were generated using the proprietary platform owned and maintained by Grzybowski Scientific Inventions, Inc., a subsidiary of Merck KGaA, Darmstadt Germany. The results are provided on an as is basis, and shall be used solely in connection with the rights afforded in the license agreement and for no other purpose.

Strategies: none selected

FGI Coeff: 0

Tunnels Coeff: 0

JSON Parameters: {}

2 Paths

5 paths found. Paths are sorted by score. Reactions are sorted in appearance order for each path.

2.1 Path 1

Score: 31.25





Figure 1: Outline of path 1

2.1.1 Appel Reaction



Substrates:

 $1. \ 1\hbox{-}(2\hbox{-bromophenyl})\hbox{-}1\hbox{-propen-}3\hbox{-ol}$

Products:

1. 1-(2-bromophenyl)-1-propen-3-yl bromide

Typical conditions: PPh3.CBr4


Protections: none

Reference: 10.1021/ja800574m and 10.1016/j.tet.2012.05.010 and

10.1016/j.tet.2004.09.021 (experimental)

Retrosynthesis ID: 9990037

2.1.2 Shapiro reaction followed by alkyl bromide addition



Substrates:

- 1. 1-(2-bromophenyl)-1-propen-3-yl bromide
- 2. Hexahydrobenzaldehyde available at Sigma-Aldrich

Products:

1. Brc1cccc1/C=C/CC=C1CCCCC1

 $\textbf{Typical conditions:} \quad 1. TsNH2NH2.2. Mes 2 Mg. LiCl. THF. heating \ then \ alkylor{1}{l} alkylor{1}{l} alkylor{2}{l} alkylo$

bromide.cooling

Protections: none

Reference: 10.1016/S0040-4039(00)75263-4 and 10.1021/ol300652k and

10.1021/jo015699l

2.1.3 Br/Li exchange

Substrates:

- $1. \ \, Brc1ccccc1/C = C/CC = C1CCCCC1$
- $2. \ \, \text{n-BuLi} \, \textbf{-} \quad \, \textit{available at Sigma-Aldrich}$

Products:

1. [Li]c1ccccc1/C=C/CC=C1CCCCC1


Typical conditions: nBuLi.or.tBuLi.THF.-78C

Protections: none

Reference: 10.1002/ejoc.201101490 and 10.1016/j.tet.2012.03.058 and 10.1016/j.tetlet.2015.01.032 and 10.1021/ja0541175 and 10.1016/j.tetlet.2016.06.123

Retrosynthesis ID: 30672

2.1.4 Addition of electrophiles to lithiated arenes/heteroarenes



Substrates:

 $1. \ [Li]c1ccccc1/C = C/CC = C1CCCCC1$

Products:

1. Oc1ccccc1/C=C/CC=C1CCCCC1

Typical conditions: B(OMe)3 then H2O2.THF

Protections: none

Reference: 10.1039/C7CC09187G (SI) and 10.1002/ejoc.201701142 and

10.1021/acscatal.6b03380 (SI,p.10) and 10.1002/chem.201702143 (SI)

Retrosynthesis ID: 10019525

2.2 Path 2

Score: 45.00





Figure 2: Outline of path 2

$\begin{tabular}{ll} \bf 2.2.1 & \bf Suzuki \ coupling \ of \ arylboronic \ pinacol \ esters \ with \ vinyl \ Bromides \\ \end{tabular}$



Substrates:

- 1. 2-Hydroxyphenylboronic acid pinacol ester available at Sigma-Aldrich
- 2. 1,2-Dibromoethylene available at Sigma-Aldrich

Products:

 $1. \ Oc1ccccc1/C=C/Br$

Typical conditions: Pd catalyst.base.solvent

Protections: none

Reference: 10.1021/cr00039a007 and $10.1007/3418_2012_32$ and 10.1021/cr0505268 and 10.1016/j.jfluchem.2016.01.018 and 10.1039/C3CS60197H

Retrosynthesis ID: 10355

2.2.2 Suzuki coupling of alkyl-9-BBNs with vinyl bromides

Substrates:

- 1. Oc1ccccc1/C=C/Br
- 2. 9-(3,3-pentamethyleneallyl)-9-borabicyclo3.3.1nonane

Products:

 $1. \ \, Oc1ccccc1/C = C/CC = C1CCCCC1$

Typical conditions: Pd catalyst.base.solvent

Protections: none

Reference: 10.1021/ja00183a048 and 10.1039/b707338k and 10.1016/j.tet.2015.05.039 and 10.1021/jo991064z and 10.1021/ol060290+ and 10.1246/bcsj.65.2863

Retrosynthesis ID: 25174

2.3 Path 3

Score: 51.25





Figure 3: Outline of path 3

2.3.1 Shapiro reaction followed by alkyl bromide addition



Substrates:

1. Hexahydrobenzaldehyde - available at Sigma-Aldrich

 $2. \ 1, 3\text{-Dibromo-1-propene} - \quad \textit{available at Sigma-Aldrich}$

Products:

 $1. \ \mathrm{Br/C}{=}\mathrm{C/CC}{=}\mathrm{C1CCCCC1}$

Protections: none

Reference: 10.1016/S0040-4039(00)75263-4 and 10.1021/ol300652k and

10.1021/jo015699l

2.3.2 Synthesis of boronic acid esters

Substrates:

1. Methoxyboronic acid pinacol ester - available at Sigma-Aldrich

 $2. \ \mathrm{Br/C}{=}\mathrm{C/CC}{=}\mathrm{C1CCCCC1}$

Products:

1. CC1(C)OB(/C=C/CC=C2CCCCC2)OC1(C)C

Typical conditions: Mg.THF.-78 C


Protections: none

Reference: 10.1039/B507900D p. 3167, 3170 and 10.1021/acs.orglett.5b01434 SI

p. S6

Retrosynthesis ID: 245548

2.3.3 Suzuki coupling of aryl bromides with alkenyl boronic acids pinacol esters



Substrates:

- $1. \ CC1(C)OB(/C{=}C/CC{=}C2CCCC2)OC1(C)C\\$
- 2. 2-Bromophenol available at Sigma-Aldrich

Products:

 $1. \ \, \text{Oc1ccccc1/C=C/CC=C1CCCCC1}$

 ${\bf Typical\ conditions:}\ {\bf Pd\ catalyst.base.solvent}$

Protections: none

Reference: 10.1021/cr00039a007 and $10.1007/3418_2012_32$ and 10.1021/cr0505268 and 10.1016/j.jfluchem.2016.01.018 and 10.1039/C3CS60197H

Retrosynthesis ID: 11057


2.4 Path 4

Score: 51.25



Figure 4: Outline of path 4

2.4.1 Bromination of vinylalanes



Substrates:

1. 4-(tert-Butyldimethylsilyloxy)-1-butyne - available at Sigma-Aldrich

Products:

1. CC(C)(C)[Si](C)(C)OCC/C=C/Br

Typical conditions: Schwartz's reagent.then.Br2


Protections: none

Reference: DOI: 10.1039/C2CC36604E (SI, page S18) AND DOI:

10.1080/00397910008087318

Retrosynthesis ID: 7405

2.4.2 Deprotection of TBS ethers



Substrates:

1. CC(C)(C)[Si](C)(C)OCC/C=C/Br

Products:

1. 4-bromobut-3-en-1-ol - available at Sigma-Aldrich


Typical conditions: TBAF.THF

Protections: none

Reference: 10.1016/j.tet.2013.01.017 and 10.1016/j.tet.2004.04.042

Retrosynthesis ID: 31010160

2.4.3 Appel Reaction





Substrates:

1. 4-bromobut-3-en-1-ol - available at Sigma-Aldrich

Products:

1. 1,4-dibromo-but-1-ene - available at Sigma-Aldrich

Typical conditions: PPh3.CBr4


Protections: none

Reference: 10.1021/ja800574m and 10.1016/j.tet.2012.05.010 and

10.1016/j.tet.2004.09.021 (experimental)

Retrosynthesis ID: 9990037

${\bf 2.4.4}\quad {\bf HWE/Wittig\ Olefination}$



Substrates:

1. Cyclohexanone - available at Sigma-Aldrich

2. 1,4-dibromo-but-1-ene - available at Sigma-Aldrich

Products:

1. Br/C=C/CC=C1CCCCC1

Typical conditions: 1.PPh3 or trialkylphosphite.2.base.aldehyde

Protections: none

Reference: 10.1002/anie.200705005 and 10.1021/ol052106a and

10.1021/jo00075a064 and 10.1021/ol3027297

2.4.5 Synthesis of vinylboronic acid from vinyl bromide

Substrates:

1. Methyl borate - available at Sigma-Aldrich

 $2.\ \mathrm{Br/C}{=}\mathrm{C/CC}{=}\mathrm{C1CCCCC1}$

Products:

1. OB(O)/C=C/CC=C1CCCCC1

 $\textbf{Typical conditions:} \ 1. nBu-Li. THF. cooling. 2. (MeO) 3B. cooling. 3. H+$

Protections: none

Reference: 10.1016/j.tet.2010.11.065 and 10.1021/ja806258v (suppl. Info)

Retrosynthesis ID: 24016

2.4.6 Suzuki coupling of aryl iodides with alkenyl boronic acids

Substrates:

1. 2-Iodophenol - available at Sigma-Aldrich

2. OB(O)/C=C/CC=C1CCCCC1

Products:

 $1. \ \, \text{Oc1ccccc1/C=C/CC=C1CCCCC1}$

 ${\bf Typical\ conditions:}\ {\bf Pd\ catalyst.base.solvent}$

Protections: none

Reference: 10.1021/cr00039a007 and $10.1007/3418_2012_32$ and 10.1021/cr0505268 and 10.1016/j.jfluchem.2016.01.018 and 10.1039/C3CS60197H

Retrosynthesis ID: 24886

2.5 Path 5

Score: 51.25





Figure 5: Outline of path 5

2.5.1 Esters reduction with LAH



Substrates:

1. C8H16O2Si

Products:

1. (e)-4-trimethylsilyl-3-buten-1-ol

Typical conditions: LiAlH4.THF.0-20 $\rm C$


Protections: none

Reference: 10.1016/j.ejmech.2019.112011 p. 5, 10 and

10.1016/j.ejmech.2020.112910 p. 3, 7

Retrosynthesis ID: 9910006

2.5.2 Appel Reaction



Substrates:

1. (e)-4-trimethylsilyl-3-buten-1-ol

Products:

1. (4-bromo-but-1-enyl)-trimethyl-silane


Typical conditions: PPh3.CBr4

Protections: none

Reference: 10.1021/ja800574m and 10.1016/j.tet.2012.05.010 and

10.1016/j.tet.2004.09.021 (experimental)

2.5.3 HWE/Wittig Olefination



Substrates:

1. Cyclohexanone - available at Sigma-Aldrich

 $2. \ \ (\text{4-bromo-but-1-enyl})\text{-trimethyl-silane}$

Products:

1. C[Si](C)(C)/C=C/CC=C1CCCCC1

 ${\bf Typical\ conditions:}\ 1. PPh3\ or\ trialkylphosphite. 2. base. aldehyde$

Protections: none

Reference: 10.1002/anie.200705005 and 10.1021/ol052106a and

10.1021/jo00075a064 and 10.1021/ol3027297

2.5.4 Iodination of Silyl Derivatives

Substrates:

- $1. \ C[Si](C)(C)/C = C/CC = C1CCCCC1$
- 2. N-Iodosuccinimide available at Sigma-Aldrich

Products:

1. I/C=C/CC=C1CCCCC1

Typical conditions: NIS. 50C. MeCN

Protections: none

Reference: DOI: 10.1016/j.tetlet.2011.02.057 or DOI: 10.1016/S0040-

4039(96)02000-X or DOI: 10.1016/S0040-4020(02)00334-4

Retrosynthesis ID: 9211

2.5.5 Suzuki coupling of arylboronic acids with vinyl iodides



Substrates:

1. 2-Hydroxyphenylboronic acid - available at Sigma-Aldrich

 $2.\ I/C{=}C/CC{=}C1CCCCC1$

Products:

 $1. \ \, \text{Oc1ccccc1/C=C/CC=C1CCCCC1}$

 ${\bf Typical\ conditions:}\ {\bf Pd\ catalyst.base.solvent}$

Protections: none

Reference: 10.1021/cr00039a007 and $10.1007/3418_2012_32$ and 10.1021/cr0505268 and 10.1016/j.jfluchem.2016.01.018 and 10.1039/C3CS60197H